

Compact patient chair ::: SAF-103-01

Features

- Solid aluminium frame construction
- Easy to assemble in short period of time: no tools needed
- Tested to support more than 400 lbs (182 kg)
- Seat height adjustment from 20" to 27" (51 to 67 cm)
- Removable and adjustable headrest
- Footrest can be adjusting from 90° to 170°
- Backrest fully adjustable from 100° to 170°
- Armrest full length and retractable
- Front, right or left patient entry
- Fabric Naugahyde upholstery
- Aseptic wipe down

Specifications

- Weight: 30 lbs (13.6 kg)
- Dimension fully extended: 71" x 21.5" (180 x 55 cm)
- Dimension folded: 8" x 21.5" x 20" (20 x 54 x 51 cm)
- Upholstery color: black
- Frame color: Sparkle silver

Options

- Carry bag (SAF-BAG1)
- Pneumatic cylinder (SAF-PNEUMATIC)
- Articulated headrest (SAF-HEADREST)
- Carry case for accessories (SAF-CASE1)

Accommodation Kit

- Includes: SAF-103(Compact Chair)
- SAF-112(Instrument Tray)
- SAF-052(Bracket & Post)




Illumination Kit

- Includes: SAF-103(Compact Chair)
- SAF-302-01(LED Light)
- SAF-052(Bracket & Post)



Supreme Kit

- Includes: SAF-103(Compact Chair)
- SAF-302-01(LED Light)
- SAF-112(Instrument Tray)
- SAF-052(Bracket & Post)



One-year warranty on manufacturing defect