

## LED light with sensor ::: SAF-302-01

Comes with Bracket & Post (SAF-052) to fix on Patient Chair (SAF-102)

### Features

- Reflection type design, using 4 high power LED's
- Daylight illumination with a light intensity up to 32,000 LUX
- Hands-free operation by sensor for on/off and control intensity
- Dimmable brightness at several levels
- Provides deep oral cavity viewing
- Lighting without glare for patient comfort
- High efficiency and low power consumption
- No cooling fan required
- Completely sealed housing & hygienic

### Specifications

- Weight: ± 1 kg (lamp only)
- Sensor switch at lamp crown: On/Off/Composite
- LED life-time: 50,000 hours (average)
- Light source: 4 x 3W LED (Cree brand)
- Power consumption: 15 Watts
- Color temperature: 6,000K
- Light intensity: 8,000-30,000Lux (adjustable)
- Illumination area: 75 x 170mm at 70cm
- 110 volts / 60Hz

### Options

- Aluminium carry case: 18x13x5.75" (46x33x15 cm) (SAF-CASE1)
- 220 volts / 50Hz
- Lamp only



**One-year warranty on manufacturing defects**