

## Dental Wireless O-Star Curing Light Cure... SAF-800-Star

### Features

Smart but powerful making a big difference  
7 modes to fully meet a variety of clinical needs  
Uniform beam, strong penetrability  
Less light loss  
Light intensity up to 3000 mW/cm<sup>2</sup>  
Color : Black

### Specifications

Size of main unit: 23mm×23mm×204mm  
Net weight of main unit: 108g  
Applied parts of the equipment: Top of main unit, Point cure lens  
Duty cycle of the equipment: 20 Sec on/20 Sec off

