

Cordless , Prophy Motor ::: SAF-PO

Features

- Wireless handpiece
- Quiet motor
- Ultra-small head
- Wireless, Stable, Convenient
- Dual-control design
- Cost-effective
- Ergonomic design

Spécifications

- Large-capacity lithium battery
- 5h continuous working
- A full charge can be used for at least 60 treatments
- Slim handpiece with anti-slip design (18mm diameter)
- Button control, Stable Speed
- 6 gear optional: 500 rpm 1000rpm 1500 rpm 2000 rpm 3000 rpm 4000 rpm
- Foot pedal control
- Non-gear shift control mode
- Speed range: 0 – 4000 rpm
- The Speed will increase according to the depth of the pedal depressed

