

Cordless LED curing light ... SAF-800-3

Features

- Light meter attached
- Built-in fan to avoid heat
- The battery can be replaced
- OLED screen, all angles visible
- Automatic memory, avoid setting frequently
- Ergonomics design, more comfortable in operation
- Three working modes with three criteria of light intensity
- Stable output of light intensity
- A full charge can be used more than 250 times continuously

Specifications

- Intensity : 850 - 1000 mw/cm²
- Wave Length : 420nm - 480nm
- Voltage : 100 -240 Vac 50/60Hz
- Weight : 210g
- Dimensions : 195 x 40 x 150 mm



Cordless LED curing light ... SAF-800-5

Features

- The battery can be replaced
- Time setting : 5, 10, 15, 20 sec
- A full charge can be used more than 500 times continuously
- Stable output of light intensity
- Low standby power consumption with forty day's standby time

Specifications

- Intensity : 850 - 1000 mw/cm²
- Wave Length : 420nm - 480nm
- Voltage : 100 -240 Vac 50/60Hz
- Weight : 145g
- Dimensions : 260 x 31 x 34 mm

