

VET series ... SAF-531-VET

Self-contained mobile veterinary dental units

Features

- Oiless, maintenance free compressor
- Built-in piezo electric ultrasonic scaler with 4 tips
- 2 automatic handpiece connections, 4 holes "Midwest"
- Control of the handpieces with variable speed disc type foot control
- Air / water syringe with indexing tip
- All tubing's are easy to clean, and instruments are autoclavable
- 4 instrument holders at working height
- 2 times 2L easy to reach bottles for the clean water system
- Double AC outlets for accessories
- Stainless steel air tank for more dry & clean air
- 4 wheels for easy movement (lockable)

Specifications

- Weight: ± 80 lbs | ± 36 kilos
- Noise level: ± 56 Db @ 1 meter
- Voltage: 110 Volts / 60 Hz
- Dimensions: 17.5" x 11" x 30.5" | 43.15 x 28 x 76.25 cm
- Air displacement compressor: 2.64 CFM
- Max pressure: 80 PSI (5.52 Bar)

Options

- LED illumination on ultrasonic scaler
- Built-in fiber optic handpiece configuration (EACH)
- Custom color (more than 75 colors)
- External quick disconnect water for auxiliaries
- Air Vacuum System with bottle
- Protection cover to protect your investment
- Custom logo on protection cover
- Turbine
- Polisher
- Additional instruments holder (each)



One-year warranty on manufacturing defects