


ULC ... SAF-531

Fully self-contained mobile dental cart

Features

- 3 automatic handpiece connections, 4 holes "Midwest"
- Control of the handpieces with variable speed disc type foot control
- Air / water syringe with indexing tip
- Deluxe high volume valve and tubing
- Deluxe low volume valve and tubing
- Solid collector 1 7/8"
- All tubing are easy to clean, and instruments are autoclavable
- 7 instrument holders at working height
- 1L easy to reach bottles for the clean water system
- 2L waste collector jar with overflow protection
- Double AC outlets for accessories "hospital grade"
- External quick disconnect water for auxiliaries such as Cavitron
- Oilless, maintenance free compressor and vacuum pump
- Stainless steel air tank for more dry & clean air
- On/off foot switch for vacuum pump
- 4 wheels for easy movement (lockable)



Specifications

- Weight: ± 95 lbs | ± 43 kilo
- Noise level: ± 56 Db @ 1 meter
- Voltage: 110 Volts / 60 Hz
- Dimensions: 17.5" x 11" x 30.5" | 43.15 x 28 x 76.25 cm
- Air displacement vacuum pump: 5.5 CFM
- 11Air displacement compressor: 2.64 CFM
- Max pressure: 80 PSI (5.52 Bar)

Options

- Amalgam separator
- Articulated arms for assistant with 2nd syringe
- ¹Automatic suction instead of pedal
- Built-in fiber optic handpiece configuration (EACH)
- Built-in piezo electric ultrasonic scaler
- Custom color



One-year warranty on manufacturing defects

¹ not compatible with articulated arms for assistant