

Expedition ... SAF-430-3
Most quiet self-contained dental unit in the market

Features

- Heavy duty "TIG" welded aluminium case
- Back pivoted handles & two 8" wheels for easy moving
- Midwest 4 hole connection for use of a high speed handpiece (wet)
- Midwest 4 hole connection for use of a low speed handpiece (dry)
- Control of the handpiece through a rheostat foot pedal
- Air / water syringe with quick change tip
- Deluxe "HVE" valve and tubing with saliva ejector tip adaptor (see option)
- All tubing are easy to clean, and instruments are autoclavable
- Removable instrument holders at adequate working height
- 1 litre bottle for the clean water system
- 2 litres waste collector jar with overflow protection
- Double AC outlets for accessories "hospital grade"
- External quick disconnect air & water for auxiliaries such as Cavitron
- Oilless, maintenance free compressor
- Stainless steel air tank for more dry air
- Oilless, maintenance free vacuum pump
- On/off foot switch for vacuum pump

Specifications

- Weight : 62.8 lbs (28.5 kg)
- Noise level : 44 Db @ 1 metre
- Voltage : 110 Volts / 60 Hz
- Dimensions : 12" x 23" x 14" (30 x 59 x 36 cm)
- Air displacement compressor : 2.3 CFM (65 L/min)
- Air displacement vacuum pump : 2.6 CFM (73.6 L/min)
- Pressure : 80 PSI (5.52 Bar)

Options

- Upper tray with 6 deluxe instruments holders (SAF-432)
- Upper tray with 2 articulated arms & 6 deluxe instruments holders (SAF-433)
- Amalgam separator
- Additional deluxe instrument holder for SAF-432 or SAF-433 (H-1075G)
- Handpiece, low speed motor, contra angle, prophylaxis ...
- Change the "HVE" valve and tubing for a Saliva Ejector valve and tubing
- Unit in 220 Volts / 50 Hz



One-year warranty on manufacturing defects

WWW.SAFARIDENTAL.COM

1 800 567.0013 | 94-D Boul. des Entreprises, Boisbriand, Qc J7G 2T3 | info@safaridental.com