

Chair mounted tray

Comes with Bracket & Post (SAF-052) to fix on Patient Chair (SAF-102)

Features

- Easy to reach & Easy setup
- Pin mount at 3 different place for smoother movement
- Two 12" arms included for precise positioning
- Adjustable in height on post (SAF-052) from 22 to 40" off ground
- Stainless steel instruments tray 13.5" x 9.75" (34 x 25 cm)(Ritter size)
- Weight : 3.6 lbs (1.7 kg)

Options

- Aluminium carry case (SAF-CASE1)
- Instrument holders (SAF-112H)



One-year warranty on manufacturing defects