

Bicarbonate jet, Handpiece ... SAF-871

Features

- Modern, compact design easy to use
- A new generation of bicarbonate jets
- Optimizing the performance of the procedures
- Container of bicarbonate easy to access
- Irrigation with pneumatic system
- Autoclavable bicarbonate jet point
- Mixes air + water+ bicarbonate at a small distance
- Transparent cover, positioned on the upper part of the pen
- Allows the quantity in the reservoir to be checked
- Resistant to oxidation Mobility in installation
- Easy to clean and disinfect
- Allows accuracy in the work and greater sensitivity
- Allows safe sterilization, avoiding cross contamination
- Ergonomics with excellent facility in handling
- Using a single foot control for activation
- Body made of light metal



Specification

- Input water pressure: 0.7bar to 2.2bar (70-220kPa)
- Input air pressure: 3.5bar to 4.5bar (350-450kPa)
- Working mode: continuous operation, Intermittent operation
- Weight: About 0.150kg
- Connector: Meets the requirements of YY / T 0514
- Operation environment:
 - Environment temperature: + 10~ +40
 - Relative humidity: 30% ~ 75% Atmospheric pressure : 50kPa ~ 106kPa