

## Dentacenter CLC with quartz top ::: SAF-530-01

Fully self-contained dental cart with powerful suction!

### Features

- 3 automatic 4 holes "Midwest" connection for handpiece
- Control of the handpiece with variable speed disc type foot control
- Air / water syringe with indexing tip
- Deluxe "HVE" valve and tubing & deluxe Saliva Ejector valve and tubing
- All tubing's are easy to clean, and instruments are autoclavable
- Line cleaner system included (chlorhexidine)
- 7 instrument holders at working height
- 2L easy to reach bottles for the clean water system
- 4L waste collector jar with overflow protection
- Double AC outlets for accessories "hospital grade"
- External quick disconnect air & water for auxiliaries such as Cavitron
- External quick disconnect water city input connection and tubing
- Oil-less, maintenance free compressor and vacuum pump
- Stainless steel air tank for more dry air
- On/off foot switch for vacuum pump
- 5 wheels for an easy movement (lockable)
- 5 instrument drawers



### Specifications

- Weight : approx 275lbs
- Noise level : 44 Db @ 1 meter
- Voltage : 110 Volts / 60 Hz
- Dimensions : 27,5" x 26" x 33"
- Air displacement vacuum pump : 5.5 CFM
- Air displacement compressor : 2.64 CFM
- Max pressure : 80 PSI (5.52 Bar)

### Options

- Built-in amalgam separator
- Tray with 2 articulated arms
- 2 Automatic suction
- Built-in fiber optic handpiece configuration (EACH)
- Built-in piezo electric ultrasonic scaler
- Additional drawers (Maximum 8 drawers total)
- Custom color (more than 75 colors)

**One year warranty on manufacturing defect**