

Halogene light ::: SAF-330

Comes with Bracket & Post (SAF-052) to fix on Patient Chair (SAF-102)

Features

- Medical halogen light
- Easy on hand positioning with flexible goose neck of 30"(76cm)
- Narrow beam of concentrated light
- 60% heat dissipation through the back of the lamp
- Double chamber reflector, fix focus
- Protective glass in case of bulb bursting

Specifications

- Weight : 3.5 lbs (1.6 kg)
- Color temperature : 3000 degree Kelvin
- Illumination area : 2.5" (6cm) diameter at 24" (61 cm)
- Intensity : 1100 candles
- Low voltage transformer: 12 volts
- Bulb : JC 12 volts 20 watts

Options

- Aluminium carry case : 18" x 13" x 5.75" (46 x 33 x 15 cm)(SAF-CASE1)
- 220 volts / 50Hz (SAF-330-220)
- Replacement bulb, 1000 hrs average (SAF-330B)
- Lamp only



One-year warranty on manufacturing defects