

Portable micro motor "E" type ... SAF-670-3

Features

About 9 hours time by full charging
Can install straight or contra angle
Rechargeable & Portable Dental Machine

Controller unit specifications

Weight : 185g
Dimensions : 67 x 132 x 30mm
Input : AC 100V ~ AC 240V, 50/60Hz
Output : DC 13V 1.2A
Charging time : About 2.5 hours

"E" Type handpiece specification

Weight : 71g
Size : 110 x Ø22mm
Speed : Max 25,000 rpm
Optional : Contra-angle or Straight nosecone

*** For sale outside Canada only**



Portable laboratory micro motor ... SAF-670-4

Features

About 9 hours time by full charging
Rechargeable & Portable Dental Machine

Controller unit specifications

Weight : 185g
Dimensions : 67 x 132 x 30mm
Input : AC 100V ~ AC 240V, 50/60Hz
Output : DC 13V 1.2A
Charging time : About 2.5 hours

Straight handpiece specification

Weight : 120g
Size : 130 x Ø24mm
Speed : Max 25,000 rpm



One-year warranty on manufacturing defect