

Assistant Stool ::: SAF-204

Features

- Adjustable abdominal support arm with comfortable shape
- Ambidextrous with its revolutionary shape
- Opening and closing mechanism ratchet
- Adjustable in height for more comfort
- Foot rest (18") adjustable in height
- Ergonomic seat cover with Naugahyde
- Pneumatic cylinder for seat height adjustment
- Backrest adjustable in height & depth
- Quick release points for easy disassemble
- Fold easily into a compact carrying bag
- Five-legs dual wheel base for stability and easy movement

Specifications

- Seat height adjustment : 16 ½" to 21" (42 to 54 cm)
- Seat size : 17 x 15" (43.18 x 38.1 cm)
- Backrest adjustment (height) : 1.5" (3.8 cm)
- Backrest adjustment (depth) : 2.5" (6.35 cm)
- Backrest size : 15 x 6" (38.1 x 15.24 cm)
- Weight : 26.4 lbs (12 kg)
- Folded dimensions : 20" x 20" x 7" (51 x 51 x 18 cm)
- Color : Black (Black US 593)
- Material : Steel and Nylon

Options

- Carry bag (SAF-BAG1)
- Without backrest (SAF-204W) (subtract)



One-year warranty on manufacturing defect