

Doctor Stool ::: SAF-203

Features

- Removable quick fit backrest
- Ergonomic seat cover with Naugahyde
- Pneumatic cylinder for seat height adjustment
- Backrest adjustable in height & depth
- Quick release points for easy disassemble
- Folds easily into a compact carrying bag
- Five-legs dual wheel base for stability and easy movement

Specifications

- Seat height adjustment : 16 ½" to 21" (42 to 54 cm)
- Seat size : 17 x 15" (43.18 x 38.1 cm)
- Backrest adjustment (height) : 1.5" (3.8 cm)
- Backrest adjustment (depth) : 2.5" (6.35 cm)
- Backrest size : 15 x 6" (38.1 x 15.24 cm)
- Weight : 15.9 lbs (7.2 kg)
- Folded dimensions : 20" x 20" x 7" (51 x 51 x 18 cm)
- Color : Black
- Material : Steel and Nylon

Options

- Carry bag (SAF-BAG2)
- Without backrest (SAF-203W) (subtract)



One-year warranty on manufacturing defects