

## **Expedition XP ... SAF-430-3XP**

**Quiet self-contained dental unit with powerful suction !**

### Features

- Heavy duty "TIG" welded aluminium case
- Back pivoted handles & two 8" wheels for easy moving
- Midwest 4 hole connection for use of a low speed handpiece (dry)
- Midwest 4 hole connection for use of a high speed handpiece (wet)
- Control of the handpieces through a rheostat foot pedal
- Air / water syringe with quick change tip
- Deluxe "HVE" valve and tubing
- Deluxe Saliva Ejector valve and tubing
- All tubing are easy to clean, and instruments are autoclavable
- Removable instrument holders at adequate working height
- 1 litre bottle for the clean water system
- 2 litres waste collector jar with overflow protection
- Double AC outlets for accessories "hospital grade"
- External quick disconnect air & water for auxiliaries such as Cavitron
- Oilless, maintenance free compressor
- Stainless steel air tank for more dry air
- Oilless, maintenance free vacuum pump
- On/off foot switch for vacuum pump

### Specifications

- Weight : 70.4 lbs (32 kg)
- Noise level : 44 Db @ 1 metre
- Voltage : 110 Volts / 60 Hz
- Dimensions : 12" x 23" x 14" (30 x 59 x 36 cm)
- Air displacement compressor : 2.3 CFM (65 L/min)
- Air displacement vacuum pump : 4.6 CFM (127.4 L/min)
- Pressure : 80 PSI (5.52 Bar)

### Options

- Only low speed handpiece controls and tubing (dry) (SAF-430-2XP)
- Amalgam separator
- Upper tray with 4 deluxe instruments holders (SAF-432)
- Upper tray with 2 articulated arms & 4 deluxe instruments holders (SAF-433)
- Additional deluxe instrument holder for SAF-432 **or** SAF-433 (H-1075G)
- Handpiece, low speed motor, contra angle, prophylaxis ...
- Unit in 220 Volts / 50 Hz



**One-year warranty on manufacturing defects**