

Hygiene unit ::: SAF-430-2

The most effective dental hygiene unit at a low price!

Features

The vacuum pump and the compressor operates independently
Heavy duty "TIG" welded aluminium case
Back pivoted handles & two 8" wheels for easy moving
Double AC outlets for accessories "hospital grade"
External quick disconnect water for auxiliaries such as Cavitron
All tubings are easy to clean, and instruments are autoclavable

Independent control of the handpiece through a rheostat foot pedal
On/off foot switch for vacuum pump
Midwest 4 holes connection for use of a low speed handpiece (dry)
Deluxe "SE" valve and tubing (see option)
Air / water syringe with quick change tip

Oilless, maintenance free compressor
Oilless, maintenance free vacuum pump
Removable instrument holders at adequate working height
1 litre bottle for the clean water system
2 litres waste collector jar with overflow protection

Specifications

Weight : 62 lbs (28.12 kg)
Noise level : 44 Db @ 1 metre
Voltage : 110 Volts / 60 Hz
Dimensions : 12" x 23" x 14" (30 x 59 x 36 cm)
Air displacement compressor : 2.3 CFM (65 L/min)
Air displacement vacuum pump : 2.6 CFM (73.6 L/min)
Pressure : 80 PSI (5.52 Bar)

Options

Upper tray with 4 deluxe instruments holders (SAF-432)
Upper tray with 2 articulated arms & 4 deluxe instruments holders (SAF-433)
Additional deluxe instrument holder for SAF-432 **or** SAF-433 (H-1075G)
Handpiece, low speed motor, contra angle, prophylaxis ...
Addition of the "HVE" valve, tubing and connector
Unit in 220 Volts / 50 Hz



One-year warranty on manufacturing defects

