

Explorer ::: SAF-500-M4

Features

- Hard case with handle for easy movement*
- Connection "Midwest" 4 holes to use low-speed handpiece (dry)*
- Connection "Midwest" 4 holes to use high speed handpiece (wet)*
- Control of the handpiece through a rheostat-type foot pedal*
- Air water syringe with quick change tip*
- Suction o"type" Venturi with saliva ejector and 7" tubing (see options)*
- All lubing are easy to clean, and instruments are autoclavable*
- 4 instrument holders included*
- 750 ml bottle for the clear water system*
- 1 litre waste reservoir*
- External quick disconnect water for auxiliaires such as Cavitror*



Specifications

- Weight : 31.1 lbs (14.06 kg)
- Dimensions : 12" x 8.5" x 22" (30.5 x 20.3 x 55.0 cm)

Options

- Change the Saliva Ejector valve and tubing for a "HVE" valve and tubing

One-year warranty on manufacturing defects