

## Piezo ultrasonic scaler with LED ::: SAF-830-7

### Features

- Automatic frequency tracking ensures that the machine always works on the best frequency and performs more steadily.
- The handpiece is detachable and can be autoclaved to the high temperature of 135°C and the pressure of 0.22MPa.
- Digitally controlled, easy operation and more efficient for scaling
- Lightweight handpiece with on/off foot control
- Titanium alloy cast handpiece assembled with excellent piezo-ceramic
- Detachable handpiece and tip are autoclavable
- The cable is made from silica gel tube, soft and durable
- Includes one of each of the following EMS-style tips: E1, G1, G2, G4, P1
- Automatic frequency tracking

### Specifications

- Power input: AC100V-240V 50/60Hz
- Main unit input: 24V 50/60Hz 1.3A (Max)
- Output power: 3-20W
- Frequency: 28kHz±3kHz
- Dimensions: 196mm \*134mm \*80mm
- Water pressure: 0.01MPa~0.5MPa



**One-year warranty on manufacturing defect**