

SAF-COCOON F

ERGONOMIC DESIGN

- Gun Grip Exposure Trigger
- Intuitive User Interface
- Easy to Clean UV Coated Case

LIGHT & COMPACT

- Weighing in at less than 5 lbs, COCOON takes radiographs twice as fast, reducing operation time and allowing users to take it anywhere.

BATTERY

- The cradle provides a convenient and immediate charging platform, while the Li-battery has a long-lasting life cycle, taking just 1 hr to be fully charged.
- In addition, the AUTO SHOW-DOWN function control power consumption to efficiently manage the battery.

ADVANCED RADIATION SAFETY

- Internally: Multiple lead wrapping ensures no radiation leakage.
- Externally: Back-scatter shield provides safety from radiation dose.
- Trigger Safety & Screen Lock: Prevents un-intended exposure

STABLE POSITIONING: AIMING LIGHT

- Rotary plate provides efficient X-ray positioning, and an Aiming light is used to accurately pinpoint the x-ray area

